

Product Datasheet

Myoglobin (Mb) Antibody HM213

Product Name	Myoglobin (Mb) Antibody
Catalog Number	HM213
Species/Host	Mouse
Tested Applications	LF
Immunogen	Native human Myoglobin
Brand	EastCoast Bio
Form/Appearance	Liquid; Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	Human Myoglobin
Matched Pairs	HM212 (capture)