

Product Datasheet

Candida albicans Antibody HM216

Product Name	Candida albicans Antibody
Catalog Number	HM216
Species/Host	Rabbit
Tested Applications	ELISA, IFA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.2-7.5
Cross Reactivity	Non-reactive with <i>C. krusei</i> , <i>Aspergillus</i> , <i>Pichia pastoris</i> , <i>Coccidioides</i> , <i>Histoplasma</i> , <i>Penicillium notatum</i> , <i>Saccharomyces cerevisiae</i>
Specificity	Candida species
Matched Pairs	HM216 (self-pairs) or JP004 (capture)
Preservatives	0.1% Sodium Azide