

Product Datasheet

Candida albicans Antibody HM216

| | |
|----------------------------|--|
| Product Name | Candida albicans Antibody |
| Catalog Number | HM216 |
| Species/Host | Rabbit |
| Tested Applications | ELISA, IFA |
| Brand | EastCoast Bio |
| Form/Appearance | Protein A Purified Monoclonal Antibody |
| Storage | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Buffer | 10 mM Phosphate Buffered Saline, pH 7.2-7.5 |
| Cross Reactivity | Non-reactive with <i>C. krusei</i> , <i>Aspergillus</i> , <i>Pichia pastoris</i> , <i>Coccidioides</i> , <i>Histoplasma</i> , <i>Penicillium notatum</i> , <i>Saccharomyces cerevisiae</i> |
| Specificity | Candida species |
| Matched Pairs | HM216 (self-pairs) or JP004 (capture) |
| Preservatives | 0.1% Sodium Azide |