

## Product Datasheet

### Procalcitonin (PCT-katacalcin) Antibody HM220

<b>Product Name</b>	Procalcitonin (PCT-katacalcin) Antibody
<b>Catalog Number</b>	HM220
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	LF
<b>Immunogen</b>	Purified recombinant human Procalcitonin
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Protein A Purified Monoclonal Antibody
<b>Storage</b>	Short Term ( $\leq 2$ weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	>95% by SDS-PAGE
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2
<b>Specificity</b>	Katacalcin (C-terminal) epitope
<b>Matched Pairs</b>	HM078 (capture)
<b>Preservatives</b>	0.09% Sodium Azide