

Product Datasheet

Luteinizing Hormone ? (LH-?) Antibody HM222

Product Name	Luteinizing Hormone α (LH- α) Antibody
Catalog Number	HM222
Species/Host	Mouse
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	2-8°C
Isotype	IgG
Clonality	Monoclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	LH 100%, β hCG < 1%, FSH < 1%, ihCG 7%, TSH < 1%
Specificity	LH α
Matched Pairs	FA313 (detection)
Preservatives	0.1% Sodium Azide