

## Product Datasheet

### Human Serum Albumin Antibody HM223

<b>Product Name</b>	Human Serum Albumin Antibody
<b>Catalog Number</b>	HM223
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA, WB
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Protein A Purified Monoclonal Antibody
<b>Storage</b>	2-8°C
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.4
<b>Cross Reactivity</b>	No reactivity with BSA, egg albumin, or human alpha fetoprotein
<b>Specificity</b>	Human Serum Albumin
<b>Preservatives</b>	0.09% Sodium Azide