

## Product Datasheet

### Parainfluenza Virus type 1 (PIV-1) Antibody HM226

<b>Product Name</b>	Parainfluenza Virus type 1 (PIV-1) Antibody
<b>Catalog Number</b>	HM226
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA, IFA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Protein A Purified Monoclonal Antibody
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Purity</b>	>90%
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2
<b>Specificity</b>	Parainfluenza Virus Type 1 & 3 fusion protein
<b>Matched Pairs</b>	HM227
<b>Preservatives</b>	0.1% Sodium Azide