

Product Datasheet

Troponin I (24-40) Antibody HM233

Product Name	Troponin I (24-40) Antibody
Catalog Number	HM233
Species/Host	Mouse
Tested Applications	LF, ELISA, TIA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified and Filtered Monoclonal Antibody
Storage	2-8°C
Isotype	IgG1, kappa
Clonality	Monoclonal
Purity	≥95%
Buffer	20 mM Sodium Phosphate, pH 7.2 with 150 mM NaCl
Specificity	Troponin I (epitope 24-40 a.a.)
Matched Pairs	HM456 (detection)