

## Product Datasheet

### Benzodiazepine (Oxazepam) Antibody HM237

<b>Product Name</b>	Benzodiazepine (Oxazepam) Antibody
<b>Catalog Number</b>	HM237
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	LF, ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Protein A Purified Monoclonal Antibody
<b>Storage</b>	Short Term ( $\leq 2$ weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Note</b>	Meets 300 ng/ml cut off for lateral flow
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>Purity</b>	>95%
<b>Buffer</b>	20 mM Phosphate Buffered Saline, pH 7.4 with 50 mM NaCl
<b>Specificity</b>	Oxazepam, Diazepam, and other Benzodiazepines
<b>Matched Pairs</b>	LA260 (capture)
<b>Preservatives</b>	0.1% Sodium Azide