

Product Datasheet

Human Chorionic Gonadotropin ? (HCG?) Antibody HM240

Product Name	Human Chorionic Gonadotropin α (HCG α) Antibody
Catalog Number	HM240
Species/Host	Mouse
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Meets 10 mIU/ml cut off
Isotype	IgG1
Clonality	Monoclonal
Purity	>95%
Buffer	20 mM Phosphate Buffered Saline, pH 7.4 with 50 mM NaCl
Cross Reactivity	No reactivity with β subunits of HCG, LH, TSH, or FSH
Specificity	HCG α
Matched Pairs	HM239 (detection)
Preservatives	0.1% Sodium Azide