

Product Datasheet

Procalcitonin (PCT) Antibody HM243

| | |
|----------------------------|---|
| Product Name | Procalcitonin (PCT) Antibody |
| Catalog Number | HM243 |
| Species/Host | Mouse |
| Tested Applications | LF, ELISA |
| Immunogen | Recombinant procalcitonin |
| Brand | EastCoast Bio |
| Form/Appearance | Protein A Purified Monoclonal Antibody |
| Storage | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Purity | >90% |
| Buffer | 10 mM Phosphate Buffered Saline, pH 8.0 with 50 mM NaCl |
| Specificity | Human procalcitonin |
| Matched Pairs | HM244 (capture) |
| Preservatives | 0.1% Sodium Azide |