

Product Datasheet

Procalcitonin (PCT) Antibody HM244

Product Name	Procalcitonin (PCT) Antibody
Catalog Number	HM244
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Recombinant procalcitonin
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	20 mM Phosphate Buffered Saline, pH 7.4 with 50 mM NaCl
Specificity	Human procalcitonin
Matched Pairs	HM243 (detection)
Preservatives	0.1% Sodium Azide