

## Product Datasheet

### Hepatitis B Surface Antigen (HBsAg) Antibody HM252

<b>Product Name</b>	Hepatitis B Surface Antigen (HBsAg) Antibody
<b>Catalog Number</b>	HM252
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	LF, ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Ion Exchange Purified Monoclonal Antibody
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Note</b>	Sensitivity of 0.5 ng/ml in Lateral Flow
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	>95%
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 8.0 with 50 mM NaCl
<b>Specificity</b>	HBsAg (including adw, ayw, and adr subtypes)
<b>Matched Pairs</b>	HM251 (detection)
<b>Preservatives</b>	0.1% Sodium Azide