

Product Datasheet

Hemoglobin (Hb) Antibody HM258

Product Name Hemoglobin (Hb) Antibody

Catalog Number HM258

Species/Host Mouse

Tested Applications LF, ELISA

Brand EastCoast Bio

Form/Appearance Ion Exchange Purified Monoclonal Antibody

Storage Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Note Sensitivity of ≤ 30 ng/ml

Isotype IgG1

Clonality Monoclonal

Purity >95%

Buffer 10 mM Phosphate Buffered Saline, pH 8.0 with 50 mM NaCl

Cross Reactivity No reactivity with hemoglobin from pig, cow, goat, sheep, horse, chicken, rabbit

Specificity Human hemoglobin

Matched Pairs HM257 (detection)

Preservatives 0.1% Sodium Azide