

## Product Datasheet

### Hemoglobin (Hb) Antibody HM258

**Product Name** Hemoglobin (Hb) Antibody

**Catalog Number** HM258

**Species/Host** Mouse

**Tested Applications** LF, ELISA

**Brand** EastCoast Bio

**Form/Appearance** Ion Exchange Purified Monoclonal Antibody

**Storage** Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

**Note** Sensitivity of  $\leq 30$  ng/ml

**Isotype** IgG1

**Clonality** Monoclonal

**Purity** >95%

**Buffer** 10 mM Phosphate Buffered Saline, pH 8.0 with 50 mM NaCl

**Cross Reactivity** No reactivity with hemoglobin from pig, cow, goat, sheep, horse, chicken, rabbit

**Specificity** Human hemoglobin

**Matched Pairs** HM257 (detection)

**Preservatives** 0.1% Sodium Azide