

Product Datasheet

Aflatoxin B1 (AFT B1) Antibody HM260

Product Name	Aflatoxin B1 (AFT B1) Antibody
Catalog Number	HM260
Species/Host	Mouse
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Meets 5 ng/ml cut off
Isotype	IgG1
Clonality	Monoclonal
Purity	>95%
Buffer	20 mM Phosphate Buffered Saline, pH 7.4 with 50 mM NaCl
Cross Reactivity	No reactivity with Aflatoxins M1 or M2
Specificity	Aflatoxin B1, G1, Q1, B2, G2
Matched Pairs	LA441 (capture)
Preservatives	0.1% Sodium Azide