

Product Datasheet

Estradiol (E2) Antibody HM275

Product Name	Estradiol (E2) Antibody
Catalog Number	HM275
Species/Host	Mouse
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Protein G Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: ≤-20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	Estradiol
Preservatives	0.03% Proclin 300