

Product Datasheet

Methcathinone Antibody HM299

Product Name	Methcathinone Antibody
Catalog Number	HM299
Species/Host	Mouse
Tested Applications	ELISA
Immunogen	Methcathinone-BSA
Brand	EastCoast Bio
Form/Appearance	DEAE Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Purity	>95% by SDS-PAGE and HPLC
Buffer	0.15 M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Specificity	Methcathinone
Preservatives	0.05% Sodium Azide