

Product Datasheet

Hepatitis B Surface Antigen (HBsAg) Antibody HM301

Product Name	Hepatitis B Surface Antigen (HBsAg) Antibody
Catalog Number	HM301
Species/Host	Mouse
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Protein G Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline pH 7.2
Specificity	HBsAg, subtypes ad & ay
Matched Pairs	FK100 (capture)
Preservatives	0.1% Sodium Azide