

Product Datasheet

Ferritin Antibody HM304

Product Name	Ferritin Antibody
Catalog Number	HM304
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Ferritin isolated from human liver
Brand	EastCoast Bio
Form/Appearance	DEAE Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2
Clonality	Monoclonal
Purity	>95% by SDS-PAGE and HPLC
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Specificity	Human Ferritin
Matched Pairs	HM303 (capture or detection)
Preservatives	0.05% Sodium Azide