

Product Datasheet

Rabies Virus Antibody HM305

Product Name Rabies Virus Antibody

Catalog Number HM305

Species/Host Mouse

Tested Applications ELISA, IFA

Brand EastCoast Bio

Form/Appearance Protein A Purified Monoclonal Antibody

Storage Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Isotype IgG2a

Clonality Monoclonal

Purity >90%

Buffer 10 mM Phosphate Buffered Saline pH 7.2

Cross Reactivity Reactive with many Rabies Virus Strains

Specificity Rabies Virus surface glycoprotein

Preservatives 0.1% Sodium Azide