

## Product Datasheet

### Rabies Virus Antibody HM305

**Product Name** Rabies Virus Antibody

**Catalog Number** HM305

**Species/Host** Mouse

**Tested Applications** ELISA, IFA

**Brand** EastCoast Bio

**Form/Appearance** Protein A Purified Monoclonal Antibody

**Storage** Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

**Isotype** IgG2a

**Clonality** Monoclonal

**Purity** >90%

**Buffer** 10 mM Phosphate Buffered Saline pH 7.2

**Cross Reactivity** Reactive with many Rabies Virus Strains

**Specificity** Rabies Virus surface glycoprotein

**Preservatives** 0.1% Sodium Azide