

Product Datasheet

Neutrophil Gelatinase-associated Lipocalin (NGAL) Antibody HM306

| | |
|----------------------------|---|
| Product Name | Neutrophil Gelatinase-associated Lipocalin (NGAL) Antibody |
| Catalog Number | HM306 |
| Species/Host | Mouse |
| Tested Applications | LF, ELISA |
| Immunogen | Purified native human NGAL |
| Brand | EastCoast Bio |
| Form/Appearance | Purified Monoclonal Antibody |
| Storage | Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Purity | $\geq 95\%$ by SDS-PAGE |
| Buffer | 10 mM Phosphate Buffered Saline, pH 7.0 |
| Specificity | NGAL |
| Matched Pairs | HM307 (detection) |
| Preservatives | 0.09% Sodium Azide |