

Product Datasheet

Neutrophil Gelatinase-associated Lipocalin (NGAL) Antibody HM306

Product Name	Neutrophil Gelatinase-associated Lipocalin (NGAL) Antibody
Catalog Number	HM306
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Purified native human NGAL
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	$\geq 95\%$ by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Specificity	NGAL
Matched Pairs	HM307 (detection)
Preservatives	0.09% Sodium Azide