

Product Datasheet

NT-proBNP (13-27) Antibody HM315

Product Name	NT-proBNP (13-27) Antibody
Catalog Number	HM315
Species/Host	Mouse
Tested Applications	LF, ELISA, Fluorescence
Immunogen	Purified recombinant human NT-proBNP
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Specificity	NT-proBNP
Matched Pairs	HM316 (detection)
Preservatives	0.09% Sodium Azide