

Product Datasheet

Calprotectin Antibody HM319

Product Name	Calprotectin Antibody
Catalog Number	HM319
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Purified Human calprotectin
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (\leq 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	Phosphate Buffered Saline, pH 7.0
Specificity	Native Human Calprotectin
Matched Pairs	HM320 (detection)
Preservatives	0.09% Sodium Azide