

Product Datasheet

Pepsinogen II (PGII) Antibody HM324

Product Name	Pepsinogen II (PGII) Antibody
Catalog Number	HM324
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Purified Human pepsinogen II
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	Phosphate Buffered Saline, pH 7.0
Specificity	PGII
Matched Pairs	HM323 (capture)
Preservatives	0.09% Sodium Azide