

Product Datasheet

Zika Virus (ZIKV-NS1) Antibody HM329

Product Name	Zika Virus (ZIKV-NS1) Antibody
Catalog Number	HM329
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Recombinant Zika Virus NS1 Protein
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1, kappa
Clonality	Monoclonal
Purity	>90% by SDS-PAGE
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	No reactivity with Dengue Virus or Chikungunya Virus
Specificity	Zika Virus
Preservatives	0.05% Sodium Azide