

Product Datasheet

Oxycodone Antibody HM343

| | |
|----------------------------|---|
| Product Name | Oxycodone Antibody |
| Catalog Number | HM343 |
| Species/Host | Mouse |
| Tested Applications | LF, ELISA |
| Immunogen | Oxycodone-BSA |
| Brand | EastCoast Bio |
| Form/Appearance | Purified Monoclonal Antibody |
| Storage | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Note | Meets 100 ng/ml in urine testing |
| Isotype | IgG1, kappa |
| Clonality | Monoclonal |
| Purity | >90% by SDS-PAGE |
| Buffer | 0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride |
| Cross Reactivity | No reactivity with codeine, hydrocodone or hydromorphone and ≤ 1% with morphine |
| Specificity | Oxycodone |
| Matched Pairs | LA266 (capture) |
| Preservatives | 0.05% Sodium Azide |