

Product Datasheet

Dengue Virus (DENV-pan NS1) Antibody HM349

Product Name	Dengue Virus (DENV-pan NS1) Antibody
Catalog Number	HM349
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Purified Dengue Virus
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2a, kappa
Clonality	Monoclonal
Purity	>90% by SDS-PAGE
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	No reactivity with Zika or Chikungunya viruses
Specificity	All Dengue Virus serotypes
Matched Pairs	HM843, HM388 or HM389 (capture)
Preservatives	0.05% Sodium Azide