

Product Datasheet

Human Immunoglobulin M (IgM) Antibody HM351

Product Name	Human Immunoglobulin M (IgM) Antibody
Catalog Number	HM351
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Human IgM
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1, kappa
Clonality	Monoclonal
Purity	>90%
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Specificity	Human IgM
Preservatives	0.05% Sodium Azide