

Product Datasheet

Hepatitis B Surface Antigen (HBsAg) Antibody HM360

Product Name	Hepatitis B Surface Antigen (HBsAg) Antibody
Catalog Number	HM360
Species/Host	Mouse
Tested Applications	LF, ELISA, WB
Immunogen	Recombinant HBsAg
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	HBsAg
Matched Pairs	HM361 (capture)
Preservatives	0.05% Proclin 300