

Product Datasheet

Hepatitis B Surface Antigen (HBsAg) Antibody HM360

| | |
|----------------------------|---|
| Product Name | Hepatitis B Surface Antigen (HBsAg) Antibody |
| Catalog Number | HM360 |
| Species/Host | Mouse |
| Tested Applications | LF, ELISA, WB |
| Immunogen | Recombinant HBsAg |
| Brand | EastCoast Bio |
| Form/Appearance | Protein A Purified Monoclonal Antibody |
| Storage | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Purity | >90% |
| Buffer | 10 mM Phosphate Buffered Saline, pH 7.4 |
| Specificity | HBsAg |
| Matched Pairs | HM361 (capture) |
| Preservatives | 0.05% Proclin 300 |