

Product Datasheet

Cystatin C Antibody HM367

Product Name	Cystatin C Antibody
Catalog Number	HM367
Species/Host	Mouse
Tested Applications	LF, ELISA, WB
Immunogen	Cystatin C native protein
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2b
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	Cystatin C
Matched Pairs	HM368 (detection)
Preservatives	0.05% Proclin 300