

Product Datasheet

Buprenorphine (BUP) Antibody HM376

| | |
|----------------------------|---|
| Product Name | Buprenorphine (BUP) Antibody |
| Catalog Number | HM376 |
| Species/Host | Mouse |
| Tested Applications | LF, ELISA |
| Immunogen | BUP-BSA |
| Brand | EastCoast Bio |
| Form/Appearance | Purified Monoclonal Antibody |
| Storage | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Isotype | IgG1, kappa |
| Clonality | Monoclonal |
| Purity | >90% |
| Buffer | 0.05M Phosphate, pH 7.15 with 0.85% Sodium Chloride |
| Specificity | BUP, norBUP and BUP glucuronide: 100% |
| Matched Pairs | LA283 or LA284 (capture) |
| Preservatives | 0.05% Sodium Azide |