

Product Datasheet

Dehydroepiandrosterone sulfate (DHEA-S) Antibody HM378

Product Name	Dehydroepiandrosterone sulfate (DHEA-S) Antibody
Catalog Number	HM378
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	DHEA-BSA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1, kappa
Clonality	Monoclonal
Purity	>90%
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% NaCl
Specificity	DHEA-3-sulfate 100%, DHEA 20%
Matched Pairs	LA285 (capture)
Preservatives	0.05% Sodium Azide