

Product Datasheet

Methadone Antibody HM384

Product Name	Methadone Antibody
Catalog Number	HM384
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Methadone-KLH
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2b, kappa
Clonality	Monoclonal
Purity	>90%
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	EDDP 10%, EMDP 5%
Specificity	Methadone
Matched Pairs	LA279 (capture)
Preservatives	0.05% Sodium Azide