

Product Datasheet

Morphine Antibody HM385

| | |
|----------------------------|---|
| Product Name | Morphine Antibody |
| Catalog Number | HM385 |
| Species/Host | Mouse |
| Tested Applications | LF, ELISA |
| Immunogen | Morphine-BSA |
| Brand | EastCoast Bio |
| Form/Appearance | Purified Monoclonal Antibody |
| Storage | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Note | Meets 300 ng/ml cutoff |
| Isotype | IgG1, kappa |
| Clonality | Monoclonal |
| Purity | >90% |
| Buffer | 0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride |
| Specificity | Morphine, Codeine, Hydrocodone 100%, Hydromorphone 90% |
| Matched Pairs | LA281 (capture) |
| Preservatives | 0.05% Sodium Azide |