

Product Datasheet

Dengue Virus (DENV-NS1 Serotype 1) Antibody HM389

Product Name	Dengue Virus (DENV-NS1 Serotype 1) Antibody
Catalog Number	HM389
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Purified Dengue Virus
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2b, kappa
Clonality	Monoclonal
Purity	>90%
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Specificity	NS1 protein of Dengue Virus type 1
Matched Pairs	HM349 (detection) and/or HM388 (capture)
Preservatives	0.05% Sodium Azide