

## Product Datasheet

### Ketamine Antibody HM390

**Product Name** Ketamine Antibody

**Catalog Number** HM390

**Species/Host** Mouse

**Tested Applications** LF, ELISA

**Brand** EastCoast Bio

**Form/Appearance** Purified Monoclonal Antibody

**Storage** Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

**Isotype** IgG2a, kappa

**Clonality** Monoclonal

**Purity** >90%

**Buffer** 0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride

**Cross Reactivity** Ketamine 100%, Norketamine 60%

**Specificity** Ketamine

**Matched Pairs** LA280 (capture)

**Preservatives** 0.05% Sodium Azide