

Product Datasheet

Mumps Virus Antibody HM401

Product Name	Mumps Virus Antibody
Catalog Number	HM401
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Purified Mumps Virus
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Specificity	Mumps Virus
Preservatives	0.02% Sodium Azide