

Product Datasheet

Measles Virus Antibody HM402

| | |
|----------------------------|---|
| Product Name | Measles Virus Antibody |
| Catalog Number | HM402 |
| Species/Host | Mouse |
| Tested Applications | LF, ELISA |
| Immunogen | Purified Measles Virus |
| Brand | EastCoast Bio |
| Form/Appearance | Purified Monoclonal Antibody |
| Storage | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Isotype | IgG |
| Clonality | Monoclonal |
| Purity | >95% by SDS-PAGE |
| Buffer | 10 mM Phosphate Buffered Saline, pH 7.0 |
| Specificity | Measles Virus nucleoprotein |
| Preservatives | 0.02% Sodium Azide |