

Product Datasheet

Helicobacter pylori Antibody HM404

Product Name	Helicobacter pylori Antibody
Catalog Number	HM404
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Purified H. pylori
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Specificity	H. pylori
Matched Pairs	HM092 (detection) or HM405 (capture)
Preservatives	0.09% Sodium Azide