

## Product Datasheet

### Human Granzyme B ELISA Kit 850.790

<b>Product Name</b>	Human Granzyme B ELISA Kit
<b>Catalog Number</b>	850.790
<b>Description</b>	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
<b>Range</b>	31.25-1000 pg/ml
<b>Sensitivity</b>	20 pg/ml
<b>Tested Applications</b>	ELISA
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Cross Reactivity</b>	No cross reactivity with Trail, TNF $\alpha$ , IFN $\gamma$ , Fas, IL-1 $\alpha$ , IL-4, IL-6, IL-8, IL-13
<b>Specificity</b>	Recognizes both natural and recombinant human Granzyme B
<b>Sample Types</b>	Serum, Plasma, Cell culture supernatant