

Product Datasheet

Influenza A Virus Antibody HM418

Product Name	Influenza A Virus Antibody
Catalog Number	HM418
Species/Host	Mouse
Tested Applications	LF
Immunogen	Inactivated Influenza A Virus
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	2-8°C
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	No reactivity with Influenza B, Adenovirus, and RSV
Specificity	Influenza A Virus Nucleoprotein
Matched Pairs	HM419 or HM423 (detection)
Preservatives	0.05% Sodium Azide