

Product Datasheet

Influenza B Virus Antibody HM424

Product Name	Influenza B Virus Antibody
Catalog Number	HM424
Species/Host	Mouse
Tested Applications	LF
Immunogen	Inactivated Influenza B Virus
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	2-8°C
Isotype	IgG2a
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	No reactivity with Influenza A, Adenovirus, and RSV
Specificity	Influenza B Virus Nucleoprotein
Matched Pairs	HM424 (self-pairs)
Preservatives	0.05% Sodium Azide