

Product Datasheet

Adenovirus Antibody HM430

| | |
|----------------------------|--|
| Product Name | Adenovirus Antibody |
| Catalog Number | HM430 |
| Species/Host | Mouse |
| Tested Applications | LF |
| Immunogen | Inactivated Adenovirus |
| Brand | EastCoast Bio |
| Form/Appearance | Protein A Purified Monoclonal Antibody |
| Storage | 2-8°C |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Purity | ≥95% by SDS-PAGE |
| Buffer | Phosphate Buffered Saline, pH 7.4 |
| Cross Reactivity | No reactivity with Influenza A, Influenza B, and RSV |
| Specificity | Hexon Protein |
| Matched Pairs | HM429 (capture) |
| Preservatives | 0.05% Sodium Azide |