

## Product Datasheet

### Human CD213a2 ELISA Kit 850.800

<b>Product Name</b>	Human CD213a2 ELISA Kit
<b>Catalog Number</b>	850.800
<b>Description</b>	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
<b>Range</b>	156-5000 pg/ml
<b>Sensitivity</b>	150 pg/ml
<b>Tested Applications</b>	ELISA
<b>Alternative Names</b>	IL-13Ra2
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Cross Reactivity</b>	No cross reactivity has been observed with TNF-a, TRAIL, IFN-g, Fas, IL-1a, IL-4, IL-6, IL-8, IL-13, Granzyme B
<b>Specificity</b>	Recognizes both natural and recombinant IL-13R $\alpha$ 2
<b>Sample Types</b>	Serum, Plasma, Cell culture supernatant