

Product Datasheet

Respiratory Syncytial Virus (RSV) Antibody HM434

| | |
|----------------------------|---|
| Product Name | Respiratory Syncytial Virus (RSV) Antibody |
| Catalog Number | HM434 |
| Species/Host | Mouse |
| Tested Applications | LF |
| Immunogen | Purified native RSV |
| Brand | EastCoast Bio |
| Form/Appearance | Protein A Purified Monoclonal Antibody |
| Storage | 2-8°C |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Purity | >95% by SDS-PAGE |
| Buffer | 10 mM Phosphate Buffered Saline, pH 7.4 |
| Cross Reactivity | No reactivity with Influenza A, Influenza B, and Adenoviruses |
| Specificity | Fusion protein of RSV strains A & B |
| Matched Pairs | HM433 (capture) |
| Preservatives | 0.05% Sodium Azide |