

Product Datasheet

Respiratory Syncytial Virus (RSV) Antibody HM438

Product Name	Respiratory Syncytial Virus (RSV) Antibody
Catalog Number	HM438
Species/Host	Mouse
Tested Applications	LF
Immunogen	Purified native RSV
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	2-8°C
Isotype	IgG2a
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	No reactivity with Influenza A, Influenza B, and Adenoviruses
Specificity	Fusion protein of RSV strains A & B
Matched Pairs	HM437 (capture)
Preservatives	0.05% Sodium Azide